

# SONOSITE LX

Optimal access - unlimited potential



FUJIFILM  
SONOSITE

secma  
medical innovation

## SEE MORE. DO MORE.

Sonosite LX is elevating clinicians' experience by making point-of-care ultrasound more inclusive and accessible in diverse environments. From the first touch, the user is welcomed to a truly exceptional heads-up workflow experience with its expansive touchscreen display that allows rotation and extension for collaboration with team members. Sonosite LX incorporates our most advanced image clarity. Brimming with educational potential, users of all levels are empowered to explore more than ever before.

### IMAGE CLARITY ELEVATED

Our largest image ever with elevated image clarity provides clinicians greater visualisation and quicker access to touch controls.



### HEADS-UP WORKFLOW

With an optimised workflow on both the user-interface and touchscreen display, Sonosite LX helps clinicians at all levels to diagnose and treat patients efficiently and with confidence.



### OPTIMAL ACCESS

As accessibility is critical in crowded bays, the system is purposely designed with a small footprint, rotating and extending display, and vertically adjustable control panel to gain optimal bedside access.

#### ADJUSTABLE HEIGHT

Adjustable height and vertically adjustable control panel to gain optimal bedside access.

#### ADAPTABLE POSITIONING

Expansive touchscreen display that rotates, tilts and extends for ease of viewing.



SCAN  
AND  
WATCH



## ENVIRONMENT INSPIRED

Sonosite LX offers accessibility with an optimal combination of our largest clinical display and an elevated heads-up workflow to improve time to diagnosis and treatment.



### DURABLE TRANSDUCER HOLDERS

Integrated holders (no more cups) for simplified clean up and ease of use.

### TRIPLE TRANSDUCER CONNECT

Three transducer ports so your most-used transducers are at the ready with the touch OF A BUTTON.

### VERSATILE STORAGE SOLUTIONS

The flexible storage system adapts to left or right side needs—gel holder, wipes holder, storage bin, and lockable drawer to secure supplies.

### CABLE MANAGEMENT SYSTEM

Cable catcher and hooks help prevent tangled transducer cables. Even at the lowest stand height, cables are less likely to drag on the floor.

### LOW-PROFILE WHEELED BASE

The low-profile wheel base helps position the system close to the bedside.



## PERIPHERAL ACCESSORIES

### BARCODE SCANNER

Save time by automatically linking patient data to the ultrasound scan. English only, where available.

### SONOSITE ECG MODULE

The 3-lead Sonosite ECG Module connects directly to Sonosite LX ultrasound system through a dedicated ECG port. It allows the system to display an ECG trace on the clinical monitor. Includes ECG lead wire set.


### NEEDLE GUIDE STARTER KITS

Improve positioning to help guide needles to an exact target on the first attempt.

### POWERPARK CORDLESS CHARGING STATION

Dock the system in a central location while the battery charges for cord-free use. Cord-free use eliminates power cords dragging on the floor, risking contamination. No power cords to clean and disinfect.

TIRED OF READING? - CONTACT THE SECMA SALES TEAM AND GET A FREE DEMO



## GUIDE NEEDLES PRECISELY


Auto Steep Needle Profiling helps in-plane needle visualisation from multiple angles with one system interaction. The clinician simply advances the needle at the angle optimal to reach the target.

### Education

No other ultrasound manufacturer is as dedicated to ensuring clinicians have access to point-of-care ultrasound education.

Sonosite LX onboard education:

- **Scan-Along Tutorials:**  
More than 150+ scan-along education tutorials for Acute Care, Anesthesia, Pain Management, OB/GYN, Procedures, and MSK
- **How-To Videos:**  
Informational videos for Imaging Basics, System Setup, Transducers, and System Use
- **Multi-tiered User Help System:**  
Help videos, getting started guides, FAQ, system control map, and information buttons for each control



## EDUCATION AT YOUR FINGERTIPS

Sonosite is dedicated to clinician education. More than 150 scan-along education tutorials are onboard for COVID-19, diagnosis, procedures, pathology, and system and transducer use.



## TRANSDUCERS - SCAN WITH CONFIDENCE

Sonosite designs, manufactures, and tests transducers in-house. The transducers aren't only durable, they're also made for superior imaging in the most demanding environments.

- IPX7 rated for immersion\*
- Drop-tested to 1 metre
- Sonosite-manufactured transducers are covered by a Sonosite 5-year standard warranty

|  |  |  |  |
|--|--|--|--|
|   |   |    |   |
| <b>C5-1</b><br>5-1 MHz Curved<br>Scan Depth: 4.7–30 cm<br><br>EXAM TYPES:<br>Abdomen, Gynaecology, Lung,<br>Nerve, Musculoskeletal, Early<br>Obstetrics, Obstetrics, Spine           | <b>C10-3</b><br>10-3 MHz Curved<br>Scan Depth: 3–18 cm<br><br>EXAM TYPES:<br>Abdomen, Gynaecology, Lung,<br>Nerve, Musculoskeletal, Spine                    | <b>P5-1</b><br>5-1 MHz Phased<br>Scan Depth: 5–35 cm<br><br>EXAM TYPES:<br>Abdomen, Cardiac, Focused<br>Cardiac, Lung, Obstetrics, Orbital,<br>Transcranial Doppler<br>(TCD) | <b>T8-3</b><br>8-3 MHz Multi<br>Scan Depth: 4–18 cm<br><br>EXAM TYPES:<br>Cardiac, Cardiac Resuscitation                             |
|   |   |    |   |
| <b>L12-3</b><br>12-3 MHz Linear<br>Scan Depth: 1.5–9 cm<br><br>EXAM TYPES:<br>Arterial, Breast, Carotid,<br>Lung, Musculoskeletal, Nerve,<br>Ophthalmic, PIV, Superficial,<br>Venous | <b>L15-4</b><br>15-4 MHz Linear<br>Scan Depth: 1.5–6 cm<br><br>EXAM TYPES:<br>Arterial, Breast, Carotid, Mus-<br>culoskeletal, Nerve, Superficial,<br>Venous | <b>L19-5</b><br>19-5 MHz Linear<br>Scan Depth: 1–6 cm<br><br>EXAM TYPES:<br>Arterial, Lung,<br>Musculoskeletal,<br>Nerve, Ophthalmic, PIV,<br>Superficial, Venous            | <b>IC10-3</b><br>10-3 MHz Curved<br>Scan Depth: 3–15 cm<br><br>EXAM TYPES:<br>Gynaecology, Early Obstetrics,<br>Obstetrics, Prostate |

### OPTIONAL NEEDLE GUIDE STARTER KIT

CIVCO Infiniti Plus™ In-Plane Ultrasound Needle Guide CIVCO AccuSITE™,  
 Out-of-Plane Ultrasound Needle Guide, CIVCO Disposable Endocavity Needle Guide

\* Protected against water immersion for  
 30 minutes at a depth of up to 1 metre.

#### DENMARK

+ 45 8710 1425  
 info@secma.dk  
 www.secma.dk

#### SWEDEN

+46 (0) 771 60 80 00  
 info@secma.se  
 www.secma.se

#### NORWAY

+47 8151 0049  
 info@secma.as  
 www.secma.as

#### NETHERLANDS

+31 (0) 347 71 50 12  
 info@secma.org  
 www.secma.org

**secma**  
 medical innovation